
# Feelo.

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**Producent:** Aicn Medical Shenzhen, Inc. Room 4C,  
Building 1, Zhongcheng Biomedical Industrial Park, No. 21  
Linhui Road, Jinsha Community, Kengzi Street, Pingshan  
District, 518118 Shenzhen, Guangdong, P.R. China

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 Authorized EU Representative  
Company name: CMC Medical Devices & Drugs S.L.  
Address: C/ Horacio Lengo N° 18, CP 29006,  
Málaga, Spain  
Web: [www.cmcmedicaldevices.com](http://www.cmcmedicaldevices.com)  
Tel: +34951214054  
Email: [Info@cmcmedicaldevices.com](mailto:Info@cmcmedicaldevices.com)

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## SET ELEMENTS

- 1 infrared thermometer
- 1 forehead cap
- 2 AAA batteries
- 1 storage bag


## Intended Use

The infrared thermometer is designed for measuring body temperature in adults and children from 3 months of age. It works by detecting thermal radiation emitted by the body—either through direct contact with the ear canal or via a non-contact measurement from the forehead (from a distance of approximately 1.5–5.0 cm). It can be used both in home settings and in medical facilities. Thanks to infrared technology, the measurement is fast and accurate. The thermometer is especially useful for children and in situations where traditional methods of temperature measurement are inconvenient or difficult. It can also be used to measure the temperature of objects.

## Main Features


- Forehead and ear temperature measurement
- Object temperature measurement option (e.g. milk or water)
- Fast results in just 1 second
- Selectable unit: °C or °F
- Bright and easy-to-read LED display
- Silent beep for discreet notifications
- Memory for the last 32 measurements
- High measurement accuracy
- Forehead measurement from a distance of 1.5–5 cm
- Fever alert (if the temperature exceeds 37.8°C, the thermometer emits about 6 short beeps to indicate elevated body temperature)
- Automatic shut-off after 60 seconds

## Unit description

1. Age division button
2. Object temperature button
3. Forehead cap
4. Probe (sensor)
5. Display screen
6. Unit body
7. Battery cover
8. On/off and measure

## Display

1. Forehead mode
2. Ear mode
3. Temperature unit
4. Baby mode (3-12 months)
5. Toddler mode (12-60 months)
6. Mode for older children and adults (>60 months)


7. Object mode
8. Memory symbol
9. Battery symbol
10. Measurement value

- **Forehead Measurement Mode:** The icon indicates that the thermometer is currently set to forehead temperature measurement mode.
- **Ear Measurement Mode:** The icon indicates that the thermometer is currently set to ear temperature measurement mode.
- **Temperature Unit:** Displays the selected temperature unit: “°C” or “°F”.
- **Measurement Value:** The result displayed after the measurement is complete.
- **Age Range:** Icons indicating the age category:
  - 3–12 months (Infant mode)
  - 12–60 months (Toddler mode)
  - Over 60 months (Child/Adult mode)
- **Battery Symbol:** Indicates that the battery level is low and replacement is required.
- **Memory Symbol:** Indicates that the device is in memory recall mode.
- **Object Mode:** Indicates that the thermometer is set to object temperature measurement mode.

## Battery Replacement

When the battery icon „“ appears on the display, the batteries must be replaced immediately. Follow the steps below:

1. Hold the thermometer in your hand as shown.
2. Firmly press the battery cover with your fingers and slide it in the direction of the arrow until it opens (refer to the illustration).
3. Remove the used batteries and insert two new AAA batteries, making sure the polarity is correct.
4. Close the battery compartment.

**Warning:**

- If you do not plan to use the thermometer for an extended period, remove the batteries.
- Dispose of used batteries in accordance with local regulations and environmental protection guidelines.
- Do not mix old and new batteries.
- Do not mix different types of batteries: alkaline, standard (carbon-zinc), or rechargeable (Ni-Cd, Ni-MH, etc.).

**Measurement Guidelines**

1. It is recommended to test the device on yourself first to become familiar with its operation.
2. To avoid inaccurate readings due to extreme temperatures, allow the thermometer to acclimate to room temperature (15°C–40°C) for at least 30 minutes before use.
3. The person being measured should also remain in room temperature conditions for at least 30 minutes to stabilize body temperature.
4. After physical exertion, rest for at least 30 minutes before taking a measurement.
5. Do not eat, drink, or exercise immediately before or during the measurement.
6. Do not take a measurement during or immediately after feeding a child.
7. If the result seems inconsistent with how the person feels or appears unusually low, repeat the measurement after 15 minutes or verify using another method for measuring core body temperature.
8. If measurement on the forehead or in the ear canal is uncomfortable, try taking the temperature behind the ear.

9. Carefully clean the sensor with an alcohol-soaked cotton swab and wait 15 minutes before using it on another person.
10. Normal body temperature varies depending on the measurement site (ear, forehead). Always measure in the same location and do not compare results from different parts of the body.
11. In the following situations, take three measurements and consider the highest value as the correct one:
  - In children under 3 years old with low immunity, where early detection of fever is critical.
  - When using the thermometer for the first time, until consistent results are achieved.
  - When the reading appears surprisingly low.
12. A basic physiological effect called vasoconstriction may occur in the early stages of fever, causing the skin to cool. As a result, the recorded temperature may appear abnormally low.
13. Human body temperature varies throughout the day and can be affected by factors such as age, gender, skin type and thickness, external temperature, and physical activity.
14. When interpreting measurement results, take natural physiological fluctuations into account. The lowest temperature usually occurs early in the morning (around 3:00–5:00 a.m.), and the highest in the afternoon (around 3:00–5:00 p.m.). In women, body temperature may rise by 0.2–0.5°C during ovulation, and is generally slightly higher than in men.
15. Before taking a measurement, make sure the skin is dry and the forehead area is free from hair.
16. When taking multiple readings, briefly move the sensor away between measurements and wait at least 5 seconds for the most accurate result.

17. Regular measurements help you better understand your normal temperature. If you notice unusual readings or changes in well-being, pay attention to symptoms and consult a doctor.
18. Always ensure the sensor is clean before each use to guarantee measurement accuracy.


## **How to Measure Temperature**

### **Turning on the Thermometer**

1. Press the „SCAN“ button to turn on the thermometer – the screen will activate and display the last recorded temperature.
2. When you hear a beep and the screen shows the symbol “---°C”, the thermometer is ready for use.
3. The device offers three measurement modes: ear / forehead / object.

### **Measuring Temperature in the Ear**

1. Remove the forehead measurement cap – the ear icon  will appear on the screen.
2. Select the appropriate age group using the button, then gently insert the probe into the ear canal and press the „SCAN“ button.
3. The thermometer will beep and display the measured temperature.
4. The next measurement can be taken after approximately 2 seconds.

**A** - For children under 12 months: Gently pull the outer part of the ear straight back.


**B** - For children over 12 months and adults: Pull the outer part of the ear back and slightly upward.

### Measuring Temperature on the Forehead

1. Attach the forehead measurement cap – the  icon will appear on the screen.
2. Select the appropriate age group using the button, then point the thermometer toward the forehead, keeping a distance of approximately 1.5 cm, and press the „SCAN” button.
3. The thermometer will beep and display the measured temperature.
4. The next measurement can be taken after about 2 seconds.


### Correct Measurement Method:

To ensure accurate results, hold the thermometer at a distance of 1.5 to 5 cm from the forehead during measurement.

### Measuring the Temperature of Objects

1. Press the object measurement button marked with the symbol „“.
2. When the „“ icon appears on the display, point the thermometer at the object and press the „SCAN“ button.
3. The thermometer will beep and display the measured temperature.
4. The next measurement can be taken after approximately 2 seconds.

### Turning Off the Device

1. Manual shutdown – press and hold the „SCAN“ button for approximately 5 seconds.
2. Automatic shutdown – the thermometer will automatically turn off after 60 seconds of inactivity.

## Changing between Celsius and Fahrenheit

- While the device is turned off, press and hold the „SCAN“ button for approximately 8 seconds until the symbol „---“ appears on the display.
- The currently selected temperature unit will be shown. Press the „SCAN“ button again to switch between °C and °F.
- If no button is pressed for about 4 seconds, the device will automatically return to the off mode.

## How to recall readings in Memory Mode

- While the device is turned off, press and hold the „SCAN“ button for approximately 4 seconds to enter memory mode – the „M“ symbol and the most recent reading will appear on the screen.
- Press the „SCAN“ button again to browse through previously saved results – from oldest to newest.
- The thermometer stores the last 32 temperature readings.

## Error Messages

| Error Message | Problem                                                 | Solution                                                                                                                                                     |
|---------------|---------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Hi</b>     | Ear mode:<br>Temperature is higher than 42.2°C (108°F). | Operate the thermometer only between the specified temperature ranges. In the event of a repeated error message, contact your retailer or Customer Services. |
| <b>Lo</b>     | Ear mode:<br>Temperature is lower than 32°C (89.6°F).   | Operate the thermometer only between the specified temperature ranges. In the event of a repeated error message, contact your retailer or Customer Services. |
| <b>Err</b>    | Sensor not welded.                                      | Contact your retailer or Customer Services.                                                                                                                  |

|                                                                                   |                     |                                             |
|-----------------------------------------------------------------------------------|---------------------|---------------------------------------------|
| <b>ErE</b>                                                                        | EEPROM is abnormal. | Contact your retailer or Customer Services. |
|  | Low power.          | Replace two new batteries.                  |

## Disposal

Do not dispose of electrical devices together with regular household waste. Use designated collection points for electronic equipment. Contact your local authorities for information on available disposal systems. Improper disposal of electronic devices — such as dumping them in landfills — may lead to the release of hazardous substances into groundwater and the food chain, posing a risk to health and the environment.

## Maintenance, Storage, and Calibration

### Maintenance

- Remove any dirt from the device housing using a soft, dry cloth.
- Clean the temperature sensor lens with a cotton swab or soft cloth moistened with alcohol. You may disinfect the device using a 70%–80% medical alcohol solution applied to a soft cloth or swab.
- If the probe tip is damaged, contact customer service.
- Do not wash the thermometer with water or detergents containing abrasives or benzene.
- Do not immerse the device in liquid.

### Storage

- Do not store or use the device in extremely high or low temperatures, high humidity, direct sunlight, near power sources, or in dusty areas, as this may lead to inaccurate measurements.
- If the thermometer will not be used for an extended period, remove the batteries.

### Calibration

- The thermometer is factory calibrated. No recalibration is required if used according to the instructions.

- If you have any doubts about measurement accuracy, contact customer service.
- Do not attempt to modify or disassemble the thermometer.

## Technical Specifications

|                                                  |                                                                                                                                                                        |
|--------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Device name</b>                               | Feelo Family Care<br>AET-R612A                                                                                                                                         |
| <b>Standard</b>                                  | The product is made under the<br>ISO80601-2-56                                                                                                                         |
| <b>Measurement mode</b>                          | Forehead mode, Ear mode and Object temperature mode. The Ear and Forehead temperature measurement modes are adjusted mode, the Object temperature mode is direct mode. |
| <b>Measurement site</b>                          | Forehead and Ear canal                                                                                                                                                 |
| <b>Measurement range (forehead and ear mode)</b> | 32.0°C-42.2°C (89.6°F-108°F)                                                                                                                                           |
| <b>Accuracy (Laboratory)</b>                     | ±0.2°C (35.0°C - 42.0°C) / ±0.4°F (95.0°F - 107.6°F),<br>±0.3°C (±0.5°F) other range                                                                                   |
| <b>Clinical repeatability</b>                    | ≤ ±0.3°C (0.5°F)                                                                                                                                                       |
| <b>Accuracy for clinical test</b>                | ≤ ±0.3°C (0.5°F)                                                                                                                                                       |
| <b>Resolution of display</b>                     | 0.1°C (0.1°F)                                                                                                                                                          |
| <b>Measurement Units</b>                         | Celsiusz (°C) i Fahrenheit (°F)                                                                                                                                        |
| <b>Measuring time</b>                            | 1s (2s later for another measurement)                                                                                                                                  |
| <b>Memories recall</b>                           | 32 measurements recall                                                                                                                                                 |
| <b>Battery</b>                                   | 2xAAA batteries                                                                                                                                                        |
| <b>Automatic shutdown function</b>               | 60s ±10s                                                                                                                                                               |
| <b>Product life</b>                              | 5 years                                                                                                                                                                |

|                                                                                                    |                                                                                                                                                                                                                                            |
|----------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Operating conditions</b>                                                                        | Ambient temperature: od 15°C do 40°C (od 59°F do 104°F)<br>Relative humidity: 10%RH - 80%RH<br>Atmospheric pressure: 70kPa - 106kPa                                                                                                        |
| <b>Storage and transportation condition</b>                                                        | Ambient temperature: od -20°C do 55°C (od -4°F do 131°F)<br>Relative humidity: 15%RH - 93%RH<br>Atmospheric pressure: 70kPa - 106kPa                                                                                                       |
| <b>Time of thermometer from highest or lowest temperature environment to operation environment</b> | 30 minutes                                                                                                                                                                                                                                 |
| <b>High temperature hint</b>                                                                       | ≥ 37.8°C (100°F)                                                                                                                                                                                                                           |
| <b>Dimensions</b>                                                                                  | 154mm × 67mm × 19mm                                                                                                                                                                                                                        |
| <b>Weight</b>                                                                                      | 57 g without batteries                                                                                                                                                                                                                     |
| <b>Note</b>                                                                                        | Not intended to be sterilized. Not for use in an Oxygen rich environment. Not intended for use in the emergency medical services environment. Not intended for use with flammable agents. Not intended for use with flammable anesthetics. |

## Safety Precautions

- Follow the maintenance guidelines outlined in this manual.
- Use the device only as intended, as described in the instructions.
- Keep the device out of reach of children.
- Children under 11 years old should not use this device without adult supervision.
- Do not use the thermometer in extreme temperatures (below 15°C / 59°F or above 40°C / 104°F) or in extreme humidity (below 10% RH or above 80% RH). Failure to comply may result in inaccurate readings.
- Do not expose the thermometer to storage conditions outside the range of -20°C / -4°F to 55°C / 131°F or 15% RH to 93% RH, as this may lead to measurement inaccuracies.
- If the temperature sensor is damaged or loose, contact customer service.
- Do not touch or blow on the infrared sensor.
- A contaminated or damaged sensor may result in inaccurate readings.
- To clean a dirty sensor, gently wipe it with a dry, soft cloth. Do not use toilet paper or paper towels, as these may scratch the sensor and affect accuracy.
- Do not immerse the thermometer in liquids. This device is not waterproof.
- Do not use a mobile phone near the operating thermometer. Portable and mobile RF communication devices may interfere with medical equipment, causing incorrect readings.
- The sensor lens is made of a durable, special material, but if damaged, stop using the thermometer immediately and contact customer service.
- Do not drop the device.
- Do not attempt to repair or modify the thermometer yourself.
- The person taking the measurement should not touch the device and the patient simultaneously.
- If the battery or protective film is accidentally

swallowed, seek medical attention immediately.

- If the thermometer has been stored in an environment different from the one in which it will be used, allow it to acclimate for about 30 minutes at room temperature before use.
- Ensure that children do not use the thermometer unsupervised, as small parts may pose a choking hazard.
- Do not service or clean the device while it is in use.

**To reduce the risk of electric shock:**

- Do not use the device while bathing.
- For more detailed instructions regarding the storage and use of the device and any included accessories, please contact the manufacturer.

## **Electromagnetic Compatibility (EMC)**

This ME EQUIPMENT or ME SYSTEM is suitable for use in both home and hospital environments.

Warning: Do not operate high-frequency surgical equipment or ME systems in MRI-shielded rooms, where electromagnetic interference levels are high.

**Warning:** Avoid placing or using the device next to or stacked with other devices, as this may result in improper operation. If such configuration is necessary, the functionality of this device and the other equipment should be observed to ensure proper operation.

Warning: Portable RF communication devices (including peripherals such as antenna cables and external antennas) should not be used closer than 30 cm (12 inches) to any part of the thermometer, including the cables specified by the manufacturer. Otherwise, device performance may be degraded.

It is recommended to maintain a safe distance from radio and microwave-emitting devices such as Wi-Fi routers, mobile phones, or microwave ovens to avoid potential interference.

To minimize the risk of electromagnetic disturbance, do not use the thermometer in close proximity to strong sources of electromagnetic radiation.

If necessary, provide a list of all cables and their maximum lengths (if applicable), transducers, and other ACCESSORIES that may be replaced by the RESPONSIBLE ORGANIZATION and could impact the ME EQUIPMENT or ME SYSTEM's compliance with Clause 7 (EMISSIONS) and Clause 8 (IMMUNITY). ACCESSORIES may be specified in general (e.g., shielded cable, load impedance) or in detail (e.g., by MANUFACTURER and DEVICE TYPE).

1. All essential instructions necessary to maintain BASIC SAFETY and ESSENTIAL PERFORMANCE regarding electromagnetic interference throughout the expected service life.
2. The manufacturer's guidelines and declaration concerning electromagnetic emissions and immunity.

**Table 1**  
**Guidance and manufacturer's declaration – electromagnetic emissions**

| Emissions test                                        | Compliance |
|-------------------------------------------------------|------------|
| RF emissions CISPR 11                                 | Group 1    |
| RF emissions CISPR 11                                 | Class B    |
| Harmonic emissions IEC 61000-3-2                      | N/A        |
| Voltage fluctuations/ flicker emissions IEC 61000-3-3 | N/A        |

**Table 2**  
**Guidance and manufacturer's declaration – electromagnetic immunity**

| Immunity test                                                                                       | IEC 60601-1-2 Test level                                                                                                            | Compliance level                                    |
|-----------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------|
| Electrostatic discharge (ESD) IEC 61000-4-2                                                         | ±8 kV contact<br>±2 kV, ±4 kV,<br>±8 kV, ±15 kV air                                                                                 | ±8 kV contact<br>±2 kV, ±4 kV,<br>±8 kV, ±15 kV air |
| Electrical fast transient/burst IEC 61000-4-4                                                       | Power supply lines: ±2 kV<br>Input/output lines: ±1 kV                                                                              | N/A                                                 |
| Surge IEC 61000-4-5                                                                                 | Line(s) to line(s): ±1 kV<br>Line(s) to earth: ±2 kV<br>100 kHz repetition frequency                                                | N/A                                                 |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | 0% 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°<br>0% 1 cycle and 70% 25/30 cycles<br>Single phase: at 0 0% 300 cycle | N/A                                                 |

|                                              |                                                                                         |                                               |
|----------------------------------------------|-----------------------------------------------------------------------------------------|-----------------------------------------------|
| Power frequency magnetic field IEC 61000-4-8 | 30 A/m<br>50Hz/60Hz                                                                     | 30 A/m<br>50Hz/60Hz                           |
| Conducted RF IEC 61000-4-6                   | 150 kHz – 80 MHz:<br>3Vrms<br>6Vrms (in ISM and amateur radio bands)<br>80% AM at 1 kHz | N/A                                           |
| Radiated RF IEC 61000-4-3                    | 10 V/m<br>80 MHz – 2,7 GHz<br>80% AM at 1 kHz                                           | 10 V/m<br>80 MHz – 2,7 GHz<br>80% AM at 1 kHz |

Note  $U_T$  is the a.c. mains voltage prior to application of the test level.

**Table 3**  
**Guidance and manufacturer's declaration – electromagnetic Immunity**

**Radiated RF  
IEC61000-4-3  
(Test specifications for  
ENCLOSURE  
PORT IMMUNITY to RF  
wireless communications  
equipment)**

| Test Frequency (MHz) | Band (MHz) | Service                                                         |
|----------------------|------------|-----------------------------------------------------------------|
| 385                  | 380-390    | TETRA 400                                                       |
| 450                  | 380-390    | GMRS 460, FRS 460                                               |
| 710<br>745<br>780    | 704-787    | Pasmo<br>LTE 13,17                                              |
| 810<br>870<br>930    | 800-960    | GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5          |
| 1720<br>1845<br>1970 | 1700-1990  | GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS |
| 2450                 | 2400-2570  | Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7            |
| 5240<br>5500<br>5785 | 5100-5800  | WLAN 802.11 a/n                                                 |

## TWO-YEAR WARRANTY

The manufacturer provides a warranty for a period of two years from the date of purchase. The warranty is valid in Poland and within the territories of official distributors. The warranty covers damages related to manufacturing defects or improper workmanship of the product.

| Modulation                          | Modulation (W) | Distance (m) | IMMUNITY TEST LEVEL (V/m) |
|-------------------------------------|----------------|--------------|---------------------------|
| Pulse modulation 18 Hz              | 1,8            | 0,3          | 27                        |
| FM $\pm$ 5 kHz deviation 1 kHz sine | 2              | 0,3          | 28                        |
| Pulse modulation 217 Hz             | 0,2            | 0,3          | 9                         |
| Pulse modulation 18 Hz              | 2              | 0,3          | 28                        |
| Pulse modulation 217 Hz             | 2              | 0,3          | 28                        |
| Pulse modulation 217 Hz             | 2              | 0,3          | 28                        |
| Pulse modulation 217 Hz             | 0,2            | 0,3          | 9                         |

Keep the proof of purchase for the duration of the warranty period. During the warranty period, you can file a claim if the product was used correctly and according to its intended purpose.

THE WARRANTY DOES NOT COVER:

Surface damages resulting from normal use.

Damages resulting from:

- a. Accidental drops or impacts, improper use, abuse or commercial use, neglect, modification of the device,
- b. Improper maintenance, lack of care for the product,
- c. Connection to an incorrect power source,
- d. Repairs or attempts to repair by a party not authorized by the manufacturer. Any attempts to open or disassemble the device, its parts, or accessories will void the warranty.

Warranty claims must be supported by evidence. The claim must be made within the warranty period. The material covering the handle is delicate, so care should be taken during use. Scratches or abrasions on the surface of the handle will be treated as surface damages resulting from normal use. Therefore, they will not be covered by the warranty. The warranty applies to components that affect the functioning of the thermometer. Proper care of the thermometer significantly increases its lifespan. Failure to follow the recommendations may result in the loss of the warranty.

**Manufacturer:** Alich Medical Shenzhen, Inc. Room 4C, Building 1, Zhongcheng Biomedical Industrial Park, No. 21 Linhui Road, Jinsha Community, Kengzi Street, Pingshan District, 518118 Shenzhen, Guangdong, P.R. China

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 Authorized EU Representative  
Company name: CMC Medical Devices & Drugs S.L.  
Address: C/ Horacio Lengo N° 18, CP 29006, Málaga, Spain  
Web: [www.cmcmedicaldevices.com](http://www.cmcmedicaldevices.com)  
Tel: +34951214054  
Email: [Info@cmcmedicaldevices.com](mailto:Info@cmcmedicaldevices.com)

The device complies with EU Regulation 2017/745.

## RoHS

RoHS stands for Restriction of Hazardous Substances Directive. It means that the product does not contain hazardous substances such as lead, mercury, cadmium, hexavalent chromium, PBB, and PBDE in amounts exceeding permissible limits.

BF type (Body Floating) is a classification used for medical devices. It indicates that a specific part of the device is intended for direct contact with the patient's body, but not with the heart.

„BF“ means that enhanced electrical insulation is provided to protect the patient from electric shock.

Medical device

Exercise caution when using the product.

The product should be protected from moisture and stored in a dry place.

## IP22

Ingress Protection (IP) rating: The first digit (2) indicates protection against solid objects larger than 12.5 mm. The second digit (2) indicates protection against vertically dripping water when the device is tilted at an angle of up to 15 degrees.

The product is recyclable and should be sorted and processed in accordance with local recycling regulations.

This device is marked in accordance with the Polish Act on Waste Electrical and Electronic Equipment and the WEEE Directive 2012/19/EU with the symbol of a crossed-out waste bin. This marking indicates that this equipment, after its service life, cannot be disposed of with other household waste. The user is required to hand it over to those conducting the collection of waste electrical and electronic equipment. Those conducting the collection, including local collection points, stores, and municipal units, create an appropriate system that enables the disposal of this equipment. Proper handling of waste electrical and electronic equipment helps to avoid harmful health and environmental consequences arising from the presence of hazardous components and improper storage and processing of such equipment.

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**Виробник:** Alicn Medical Shenzhen, Inc. Room 4C, Building 1, Zhongcheng Biomedical Industrial Park, No. 21 Linhui Road, Jinsha Community, Kengzi Street, Pingshan District, 518118 Shenzhen, Guangdong, P.R. China

**Імпортер:** SHOP DENT Sp. z o.o.  
ul. Dziupli 12, 02-454 Warsaw, Poland  
mail: reklamacje@shop-dent.pl  
www.feelo.pl  
Модель: AET-R612A

 Authorized EU Representative  
Company name: CMC Medical Devices & Drugs S.L.  
Address: C/ Horacio Lengo N° 18, CP 29006,  
Málaga, Spain  
Web: www.cmcmedicaldevices.com  
Tel: +34951214054  
Email: Info@cmcmedicaldevices.com

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